

# EnviroScreen

## Biological Indicators

- Do you use VHP (vaporized Hydrogen peroxide) to disinfect your unit?
- Biological indicators are the gold standard in sterilisation confirmation.
- *Geobacillus thermophilus* impregnated discs supplied to you for sterilization confirmation.
- Discs are incubated for 24 h (rapid) or 5 days (full confirmation).
- We give you a validation success report, giving you confidence that the VHP was successful and the unit/cage etc is safe for use.
- We can also validate autoclaves.



Contact us at: [info@sdlab.co.uk](mailto:info@sdlab.co.uk)  
Tel: +44(0) 1483 268 300

