

## Mouse Bead Testing Profiles

Agent	SDL Mouse Quick	SDL Mouse Quarterly	SDL Mouse Annual
MHV	✓	✓	✓
MVM	✓	✓	✓
MPV	✓	✓	✓
NS-1	✓	✓	✓
Sendai		✓	✓
PVM		✓	✓
REO III		✓	✓
TMEV	✓	✓	✓
EDIM	✓	✓	✓
MNV	✓	✓	✓
M. pulmonis		✓	✓
C. piliforme		✓	✓
MAD1			✓
MAD2			✓
LCMV			✓
MCMV			✓
Ectromelia			✓
Hantaan			✓
Polyoma			✓
K-virus			✓
CARB			✓
ECUN			✓
MTLV			✓
LDV			✓
PHV			✓

## Rat Bead Testing Profiles

Agent	SDL Rat Quick	SDL Rat Quarterly	SDL Rat Annual
SDAV	✓	✓	✓
KRV	✓	✓	✓
H1	✓	✓	✓
RPV	✓	✓	✓
RMV	✓	✓	✓
NS-1	✓	✓	✓
P. carinii	✓	✓	✓
RTV	✓	✓	✓
Sendai (PIV-1)		✓	✓
PVM	✓	✓	✓
REO III		✓	✓
Hantaan		✓	✓
C. piliforme		✓	✓
M. pulmonis	✓	✓	✓
LCMV			✓
MAD1			✓
MAD2			✓
CAR Bacillus			✓
SHFV (PIV-3)			✓
E. cuniculi			✓
IDIR			✓

## Guinea Pig and Rabbit Bead Testing Profiles

Agent	SDL Guinea Pig	SDL Rabbit
Sendai (PIV-1)	✓	✓
SV-5 (PIV-2)	✓	✓
SHFV (PIV-3)	✓	
PVM	✓	
REO III	✓	✓
GPCMV	✓	
M. pulmonis	✓	
MAD1 + MAD2	✓	
LCMV	✓	✓
E. cuniculi		✓
Rabbit Rotavirus		✓
CAR Bacillus		✓
C. piliforme		✓
Toxoplasma		✓
RHDV (ELISA)		✓